

Protocol for Preparation of Coomassie G-250-stained 1D Gel Bands for In-gel Trypsin Digestion

I. Washing/Destaining

- a. The SDS-PAGE gel should only be stained with high quality Coomassie Brilliant-Blue G-250. Use only MilliQ water for destaining and reagent preparation.
- b. Wash 1.5 mL Eppendorf microcentrifuge tubes with 50% HPLC grade acetonitrile:MilliQ water.
- c. Excise gel band(s), cut into 1mm pieces and place into a washed tube(s).
- d. Equilibrate gel pieces with 500 μ L 100 mM ammonium bicarbonate (aq) buffer for 20 minutes at room temperature with gentle agitation.
- e. Wash gel pieces with 500 μ L wash solution (50 mM ammonium bicarbonate, 50% acetonitrile) until Coomassie is removed (3 - 4 x 15 min, with gentle agitation).
- f. Rinse gel pieces briefly with 500 μ L 100% HPLC grade acetonitrile. Discard acetonitrile, and allow gel pieces to air-dry.

II. Reduction/Alkylation

- a. Reduce the in-gel protein with dithiothreitol (150 μ L 10mM DTT/ 100mM ammonium bicarbonate) for 30 minutes at 56° C.
- b. Cool the sample to room temperature, remove and discard DTT solution.
- c. Alkylate the in-gel protein with iodoacetamide (100 μ L 50 mM IAA/100 mM ammonium bicarbonate) in the dark at room temperature for 30 minutes. Remove and discard IAA solution.
- d. Wash sample at room temperature for 15 minutes with 500 μ L of wash solution.
- e. Rinse gel pieces briefly with 500 μ L 100% HPLC grade acetonitrile. Discard acetonitrile and allow to air-dry.

III. Digestion/Extraction

- a. Rehydrate gel pieces for 60 minutes at 4°C (on ice) in 30 μ L of 0.02 μ g/ μ L trypsin (Promega, modified porcine, TPCK-treated, Cat#: V5111) in 50mM

ammonium bicarbonate. (Be sure that enough volume is added to ensure complete rehydration of gel pieces)

- b. Remove trypsin solution and replace with just enough 50mM ammonium bicarbonate to cover gel pieces. Incubate at 37°C for 16-18 hours (overnight).
- c. Add 7 μ L of 1% formic acid to stop the action of trypsin. Spin down the sample and transfer supernatant (containing tryptic peptides) to a clean 0.5 mL Eppendorf microcentrifuge tube and keep supernatant on ice.
- d. Rinse the gel pieces with a small volume of 0.1% formic acid and transfer solution to tube from step IIIc.
- e. Store the samples at -20°C until ready for analysis.

Notes

- Always use non-latex gloves when handling samples, keratin and latex proteins are potential sources of contamination.
- Never re-use any solutions, abundant proteins will partially leach out and contaminate subsequent samples.

References

1. Havlis J, Thomas H, Sebela M, Shevchenko A: Fast-response proteomics by accelerated in-gel digestion of proteins. *Analytical Chemistry* **75**: 1300-1306 (2003).
2. Hellman U, Wernstedt C, Góñez J, Heldin C-H: Improvement of an "in-gel" digestion procedure for the micropreparation of internal protein fragments for amino acid sequencing. *Analytical Biochemistry* **224**: 451-455 (1995).
3. Jiménez CR, Huang L, Qiu Y, Burlingame AL: In-gel digestion of proteins for MALDI-MS fingerprint mapping. In: Coligan JE (ed) *Current Protocols in Protein Science*, pp. 16.4.1-16.4.5. John Wiley & Sons, Inc., Brooklyn, N.Y. (1998).
4. Rosenfeld J, Capdevielle J, Guillemot JC, Ferrara P: In-gel digestion of proteins for internal sequence analysis after one- or two-dimensional gel electrophoresis. *Analytical Biochemistry* **203**: 173-179 (1992).
5. Shevchenko A, Wilm M, Vorm O, Mann M: Mass spectrometric sequencing of proteins from silver-stained polyacrylamide gels. *Analytical Chemistry* **68**: 850-858 (1996).