

# Protocol for Preparation of Silver-Stained 1D Protein Bands for In-Gel Tryptic Digestion

## I. Washing/Destaining

- a. Wash 0.5 mL Eppendorf microcentrifuge tubes with 50% HPLC grade acetonitrile, 0.1% trifluoroacetic acid (TFA).
- b. Using a freshly stained gel, cut the protein bands at the margin of the visible bands using a new razor blade. Wash the razor blade after cutting each band. Cut gel band into 1mm cubes and place gel pieces in pre-washed 0.5 mL microcentrifuge tube.
- c. Add 200  $\mu$ L MilliQ water and shake for 1 hour to remove the acetic acid. Remove and discard supernatant.
- d. Destain by adding 25  $\mu$ L of 15mM potassium ferricyanide and 25  $\mu$ L of 50mM sodium thiosulfate. Shake gently until the band pieces go clear (**i.e. until all the silver is removed**). Remove and discard the supernatant. The gel pieces will now appear yellow in color.
- e. Add 200  $\mu$ L of 200mM ammonium bicarbonate, pH 8. Shake for 15 minutes. Gel pieces should be colorless. Remove and discard supernatant.
- f. Add 200  $\mu$ L of 50:50 (v/v) 200 mM ammonium bicarbonate: acetonitrile. Shake gently for 15 minutes. Remove and discard supernatant.
- g. Dehydrate gel piece by adding 200  $\mu$ L HPLC grade acetonitrile. Remove immediately and discard supernatant. Allow residual acetonitrile to evaporate completely.

## II. Reduction/Alkylation

- a. Reduce the in-gel protein with dithiothreitol (150  $\mu$ L 10mM DTT in 100mM ammonium bicarbonate) for 30 minutes at 56°C.
- b. Allow sample to cool to room temperature, remove and discard DTT solution.
- c. Alkylate the in-gel protein with iodoacetamide (100  $\mu$ L 50mM IAA in 100mM ammonium bicarbonate) in the dark at room temperature for 30 minutes. Remove and discard IAA solution.
- d. Add 200  $\mu$ L of 50:50 (v/v) 200 mM ammonium bicarbonate: acetonitrile. Shake gently for 15 minutes. Remove and discard supernatant.

- e. Dehydrate gel piece by adding 200  $\mu\text{L}$  HPLC grade acetonitrile. Remove immediately and discard supernatant. Allow residual acetonitrile to evaporate completely.

### III. Digestion/Extraction

- a. Rehydrate each gel piece at 4°C in 10-20  $\mu\text{L}$  0.02  $\mu\text{g}/\mu\text{L}$  of Promega Sequencing Grade Modified Trypsin prepared in 50mM ammonium bicarbonate. If additional solvent is required to cover gel piece, use 50 mM ammonium bicarbonate (NO Trypsin).
- b. Close the microcentrifuge tubes tightly and incubate overnight at 37°C.
- c. Add 7  $\mu\text{L}$  of 1% formic acid to stop the action of trypsin. Spin down the sample and transfer supernatant (containing tryptic peptides) to a clean 0.5 mL Eppendorf microcentrifuge tube and keep supernatant on ice.
- d. Rinse the gel pieces with a small volume of 0.1% formic acid and transfer solution to tube from step IIIc.
- e. Store the samples at -20°C until ready for analysis.

### Notes

- Digestion protocol may be conveniently stopped when the water is added before shaking (step Ic). Freeze the bands/spots solution at -20°C. When ready to continue digestion, thaw bands/spots and shake for 1 hour.
- Always use non-latex gloves when handling samples, keratin and latex proteins are potential sources of contamination.
- Never re-use any solutions, abundant proteins will partially leach out and contaminate subsequent samples.

### References

1. Gharahdaghi F, Weinberg CR, Meagher DA, Imai BS, Mische SM: Mass spectrometric identification of proteins from silver-stained polyacrylamide gel: A method for the removal of silver ions to enhance sensitivity. *Electrophoresis* **20**: 601-605 (1999).
2. Havlis J, Thomas H, Sebela M, Shevchenko A: Fast-response proteomics by accelerated in-gel digestion of proteins. *Analytical Chemistry* **75**: 1300-1306 (2003).

3. Hellman U, Wernstedt C, Góñez J, Heldin C-H: Improvement of an "in-gel" digestion procedure for the micropreparation of internal protein fragments for amino acid sequencing. *Analytical Biochemistry* **224**: 451-455 (1995).
4. Jiménez CR, Huang L, Qiu Y, Burlingame AL: In-gel digestion of proteins for MALDI-MS fingerprint mapping. In: Coligan JE (ed) *Current Protocols in Protein Science*, pp. 16.4.1-16.4.5. John Wiley & Sons, Inc., Brooklyn, N.Y. (1998).
5. Rosenfeld J, Capdevielle J, Guillemot JC, Ferrara P: In-gel digestion of proteins for internal sequence analysis after one- or two-dimensional gel electrophoresis. *Analytical Biochemistry* **203**: 173-179 (1992).
6. Shevchenko A, Wilm M, Vorm O, Mann M: Mass spectrometric sequencing of proteins from silver-stained polyacrylamide gels. *Analytical Chemistry* **68**: 850-858 (1996).